

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

PRODUCT NAME: VACU-TEC

Code	10200		
Lot No.	3209	2024-05	2023-07
Reference IO-09-88			

ON LINE QUALITY CONTROLS (model CDP not enclosed)**QC Passed**

- sodium citrate solution concentration
- sodium citrate volume
- vacuum level (*)
- correct placement of the cap
- no extra body inside the product

☒
☒
☒
☒
☒
☒

(*) Only for VACU-CODE and VACU-TEC tubes

PACKAGING QUALITY CONTROL (model CDP not enclosed)**QC Passed**

- control of the number of the product for each pack

☒**QUALITY CONTROLS ON FINISHED PRODUCT** (model CDP not enclosed)**QC Passed**

- pack integrity
- data printed on the product label

☒
☒

QC Passed

N/A

Irradiation treatment Certificate n° date

☐☒**QC Passed**

N/A

**Gamma rays treatment process in compliance with the specifications
of the supply Agreement**

☐☒

Document number:

ODV N.A.

*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED

RESPONSIBLE HEAD, Quality Control Reagents.....

DATE 09/08/23