

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020

DIESSE

Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

QUALITY CONTROL CERTIFICATE

KIT NAME TOXOPLASMA IgG - CHORUS

Code 81040

Lot No. 936

2024-05

2023-05

REAGENTS

Label	Description	Lot.	Exp.	NO
D2271	TOXOPLASMA IgG	936	2024-05	6X6
B2272	CALIBRATORE TOXOPLASMA IgG	936	2024-05	1
C2271	CONTROL SERUM TOXOPLASMA IgG	936	2024-05	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	339	≥ 100
Control Serum (IU/mL)	23,6	19,4-38,1
Negative Control (IU/mL)	<5,0	<5


CALIBRATOR

0 8 2 9 3 6 6 3 2 6 4 7 3 5 5 6 0 0

Label checked with barcode reader

YES ☒

Operator Signature



CONTROL SERUM

0 8 2 9 3 6 6 2 0 3 7 0 0 0 0 0 0 1

Label checked with barcode reader

YES ☒

Operator Signature



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *Alciandri*

DATE 09/08/2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.