

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

PRODUCT NAME: VES TEC

Code 10214

Lot No. 3823



2024-11



2023-09

Reference IO-09-089

ON LINE QUALITY CONTROLS (model CDP not enclosed)

QC Passed

- sodium citrate solution concentration
- sodium citrate volume
- vacuum level (*)
- correct placement of the cap
- no extra body inside the product



(*) Only for VACU-CODE and VACU-TEC tubes



PACKAGING QUALITY CONTROL (model CDP not enclosed)

QC Passed

- control of the number of the product for each pack



QUALITY CONTROLS ON FINISHED PRODUCT (model CDP not enclosed)

QC Passed

- pack integrity
- data printed on the product label



QC Passed

N/A

Irradiation treatment Certificate n°



QC Passed

N/A

Gamma rays treatment process in compliance with the specifications of the supply Agreement



Document number:

ODV N.A.

*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED

RESPONSIBLE HEAD, Quality Control Reagents.....

DATE25/08/23