

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020

DIESSE

Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

QUALITY CONTROL CERTIFICATE

KIT NAME CALPROTECTIN - CHORUS

Code 86600

Lot No. 970

2024-09

2023-09

REAGENTS

Label	Description	Lot	Exp.	NO
D9500	CALPROTECTIN	970	2024-09	2 X 6
B9500	CALIBRATORE CALPROTECTIN	970	2025-01	1
C9500	CONTROL CALPROTECTIN	970	2025-01	1
G9502	SAMPLE DILUENT	970	2025-09	1

CONTROLS


	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	380	≥ 100
Control (mg/kg)	166,0	95,5-241,7

CALIBRATOR

9	0	0	9	7	0	6	2	2	1	2	6	3	5	3	8	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---


Label checked with barcode reader YES ☒ Operator Signature *MB*



CONTROL SERUM

9	0	0	9	7	0	6	2	0	3	7	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES ☒ Operator Signature *MB*



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents. *[Signature]*

DATE *27/03/2023*

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.