

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

QUALITY CONTROL CERTIFICATE

KIT NAME CHORUS suPAR EXTENDED

Code 81414

Lot No. 970

2024-09

2023-09

REAGENTS

Label	Description	Lot.	Exp	No
D4530	suPAR EXTENDED	970	2024-09	6X6
B4530	Calibratore suPAR EXTENDED	970	2024-09	1
C4530	Control suPAR EXTENDED	970	2024-09	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	700	≥ 100
Control (ng/ml)	8,5	6,8-12,5

CALIBRATOR

1	7	3	9	7	0	6	1	1	6	2	7	5	4	4	3	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES ☒

Operator Signature *HB*




CONTROL SERUM

1	7	3	9	7	0	6	1	4	2	3	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES ☒

Operator Signature *HB*



**This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.**

APPROVED RESPONSIBLE HEAD, Quality Control Reagents *Am*

DATE 17/10/2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.