

DIESSE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale  
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

### QUALITY CONTROL CERTIFICATE

KIT NAME INFLUENZA B IgA -CHORUS

Code 81195

Lot No. 970



2024-09



2023-10

### REAGENTS

| Label | Description                   | Lot. | Exp.    | NO  |
|-------|-------------------------------|------|---------|-----|
| D7541 | INFLUENZA B IgA               | 970  | 2024-09 | 6X6 |
| B7541 | CALIBRATORE INFLUENZA B IgA   | 970  | 2024-10 | 1   |
| C7541 | CONTROL SERUM INFLUENZA B IgA | 970  | 2024-10 | 1   |

### CONTROLS

|  | Found results | Expected results |
|--|---------------|------------------|
|--|---------------|------------------|

|                        |     |            |
|------------------------|-----|------------|
| Calibrator (mOD)       | 340 | $\geq 100$ |
| Control Serum (Ndx)    | 2,9 | 2,3-4,2    |
| Negative Control (Ndx) | 0,1 | $\leq 0,7$ |

### CALIBRATOR


|   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| 0 | 7 | 1 | 9 | 7 | 0 | 6 | 0 | 0 | 0 | 0 | 0 | 3 | 1 | 0 | 0 | 0 | 0 |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|

Label checked with barcode reader

YES



Operator Signature *FB*



### CONTROL SERUM

|   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| 0 | 7 | 1 | 9 | 7 | 0 | 6 | 2 | 3 | 4 | 2 | 0 | 0 | 0 | 0 | 0 | 0 | 1 |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|

Label checked with barcode reader

YES



Operator Signature *FB*



*This batch has been manufactured and checked according to quality requirements.  
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents *FB*

DATE 03/11/2023

### WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.