

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020

DIESSE

Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

QUALITY CONTROL CERTIFICATE

KIT NAME ANA-Screen - CHORUS

Code 86014

Lot No. 974



2024-09



2023-09

REAGENTS

Label	Description	Lot.	Exp.	NO
D3141	ANA-Screen	974	2024-09	6X6
B3141	CALIBRATORE ANA-Screen	974	2024-09	1
C3141	CONTROL SERUM ANA-Screen	974	2024-09	1

CONTROLS


	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	260	≥ 100
Control Serum (Ndx)	3,6	2,6-4,9

CALIBRATOR

1 1 5 9 7 4 6 0 0 0 0 2 7 0 0 0 0

Label checked with barcode reader YES ☒ Operator Signature RB



CONTROL SERUM

1 1 5 9 7 4 6 2 6 4 9 0 0 0 0 0 0 1

Label checked with barcode reader YES ☒ Operator Signature RB



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents [Signature]

DATE 10/11/2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.