

DIESSE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale  
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

### QUALITY CONTROL CERTIFICATE

**KIT NAME** BETA 2-GLYCOPROTEIN-M - CHORUS

Code 86052

Lot No. 988



2024-11



2023-11

### REAGENTS

Label	Description	Lot.	Exp.	NO
D3371	BETA 2-GLYCOPROTEIN-M	988	2024-11	6X6
B3371	CALIBRATORE BETA 2-GLYCOPROTEIN-M	988	2024-11	1
C3371	CONTROL SERUM BETA 2-GLYCOPROTEIN-M	988	2024-11	1

### CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	430	$\geq 100$
Control Serum (AU/mL)	98,9	50,0-157,1

### CALIBRATOR

2	0	5	9	8	8	6	3	2	4	3	9	2	9	6	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES



Operator Signature

FB



### CONTROL SERUM

2	0	5	9	8	8	6	2	4	4	5	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES



Operator Signature

FB



*This batch has been manufactured and checked according to quality requirements.  
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

FB

DATE 15/12/2023

### WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.