

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020

DIESSE

Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

QUALITY CONTROL CERTIFICATE

KIT NAME ANTI-CCP - CHORUS

Code 86094

Lot No. 989



2024-11



2023-12

REAGENTS

Label	Description	Lot.	Exp.	NO
D8671	ANTI-CCP	989	2024-11	6x6
B8671	CALIBRATORE ANTI-CCP	989	2024-11	1
C8671	CONTROL SERUM ANTI-CCP	989	2024-11	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	340	≥ 100
Control Serum (AU/mL)	68,7	46,1-90,4

CALIBRATOR


2 5 6 9 8 9 6 3 2 8 5 4 2 7 7 2 0 0

Label checked with barcode reader

YES ☒

Operator Signature

FB



CONTROL SERUM


2 5 6 9 8 9 6 2 6 4 9 0 0 0 0 0 0 1

Label checked with barcode reader

YES ☒

Operator Signature

FB



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents.....*FB*

DATE 09/01/2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.