

QUALITY CONTROL CERTIFICATE

KIT NAME CHLAMYDOPHILA PNEUMONIAE IgA - CHORUS

Code 81252

Lot No. 988



2024-11



2023-11

REAGENTS

| Label | Description                         | Lot. | Exp.    | NO  |
|-------|-------------------------------------|------|---------|-----|
| D8610 | Chlamydomphila pneumoniae IgA       | 988  | 2024-11 | 6X6 |
| B8610 | CALIBRATORE Chlamydomphila p. IgA   | 988  | 2024-11 | 1   |
| C8610 | CONTROL SERUM Chlamydomphila p. IgA | 988  | 2024-11 | 1   |

CONTROLS

|  | Found results | Expected results |
|--|---------------|------------------|
|--|---------------|------------------|

|                          |       |             |
|--------------------------|-------|-------------|
| Calibrator (mOD)         | 580   | $\geq 100$  |
| Control Serum (AU/mL)    | 38,7  | 22,2-64,0   |
| Negative Control (AU/mL) | <10,0 | $\leq 10,0$ |

CALIBRATOR

|   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| 1 | 9 | 9 | 9 | 8 | 8 | 6 | 8 | 1 | 2 | 2 | 1 | 5 | 3 | 3 | 1 | 0 | 0 |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|

Label checked with barcode reader

YES



Operator Signature

MB



CONTROL SERUM

|   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| 1 | 9 | 9 | 9 | 8 | 8 | 6 | 1 | 3 | 2 | 5 | 0 | 0 | 0 | 0 | 0 | 0 | 1 |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|


Label checked with barcode reader

YES



Operator Signature

MB



This batch has been manufactured and checked according to quality requirements.  
Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

*[Signature]*

DATE

11/01/2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.