



QUALITY CONTROL CERTIFICATE

KIT NAME HERPES SIMPLEX 1+2 IgG -CHORUS

Code 81020

Lot No. 003



2025-04



2024-01

REAGENTS

Label	Description	Lot.	Exp.	NO
D2212	Herpes simplex 1+2 IgG	003	2025-04	6X6
B2213	CALIBRATORE Herpes simplex 1+2 IgG	003	2025-04	1
C2212	CONTROL SERUM Herpes simplex 1+2 IgG	003	2025-04	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	378	≥ 100
Control Serum (INDEX)	3,1	2,1-3,8
Negative Control (INDEX)	0,0	$< 0,7$


CALIBRATOR

0 0 4 0 0 3 7 0 0 0 0 0 3 4 0 0 0 0

Label checked with barcode reader

YES ☒

Operator Signature *88*




CONTROL SERUM

0 0 4 0 0 3 7 2 1 3 8 0 0 0 0 0 0 1

Label checked with barcode reader

YES ☒

Operator Signature *88*



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *Abene*

DATE *21.03.2024*

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.