



QUALITY CONTROL CERTIFICATE

KIT NAME MEASLES IgG CHORUS

Code 81072

Lot No. 024

2025-07

2024-04

REAGENTS

Label	Description	Lot.	Exp.	NO
D2431	MEASLES IgG	024	2025-07	6X6
B2432	CALIBRATORE MEASLES IgG	024	2025-07	1
C2431	CONTROL SERUM MEASLES IgG	024	2025-07	1

CONTROLS


	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	380	≥ 100
Control Serum (Index)	3,3	2,3-4,4
Negative Control (Index)	0,0	$< 0,7$

CALIBRATOR

0 0 9 0 2 4 7 0 0 0 0 0 3 0 0 0 0 0

Label checked with barcode reader YES ☒ Operator Signature *HB*



CONTROL SERUM

0 0 9 0 2 4 7 2 3 4 4 0 0 0 0 0 0 1

Label checked with barcode reader YES ☒ Operator Signature *HB*



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents: *[Signature]*

DATE *18.06.2024*

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.